

# maxxeus SPINE





Cervical Spacer



Cervical Spacer  
with Cancellous



Cortical Disc



Cloward Dowel




Cancellous Block

Item Code	Freeze Dried	Frozen	Allograft Description	Size
1554	-12		Cervical Spacer, Lordotic	5mm
1555	-12		Cervical Spacer, Lordotic	6mm
1556	-12		Cervical Spacer, Lordotic	7mm
1557	-12		Cervical Spacer, Lordotic	8mm
1558	-12		Cervical Spacer, Lordotic	9mm
1559	-12		Cervical Spacer, Lordotic	10mm
1566	-12		Cervical Spacer, Lordotic, w/Cancellous	5mm
1567	-12		Cervical Spacer, Lordotic, w/Cancellous	6mm
1568	-12		Cervical Spacer, Lordotic, w/Cancellous	7mm
1569	-12		Cervical Spacer, Lordotic, w/Cancellous	8mm
1570	-12		Cervical Spacer, Lordotic, w/Cancellous	9mm
1571	-12		Cervical Spacer, Lordotic, w/Cancellous	10mm
1230	-12		Cortical Disc (1.2mm x 5-10mm in height)	Specify Size
1030	-12	-14	Cloward Dowel	Specify Size
1022	-12	-14	Cloward Dowel, Patellar	Specify Size
1000	-12		Cancellous Block, Cervical	Specify Size
2014-40			Maxxeus DBM Putty	0.5 cc
2015-40			Maxxeus DBM Putty	1.0 cc
2016-40			Maxxeus DBM Putty	2.5 cc
2017-40			Maxxeus DBM Putty	5.0 cc
2018-40			Maxxeus DBM Putty	10.0 cc

Cervical Spacer  
Foot Print:  
11mm x 14mm  
8° of Lordosis

Cancellous is  
Press Fit within  
the cortical  
structure



Maxxeus DBM Putty

**maxxeus**  
SPINE

Center for Tissue, Innovation and Research

2900 College Drive

Kettering, Ohio 45420-2972

1-855-MAXXEUS (629-9387)