


# COPAL<sup>®</sup> knee moulds

## Sizes and dimensions

	Femoral component		Tibial component				
Size	Media Lateral (ML)	Anterior Posterior (AP)	Media Lateral (ML)	Anterior Posterior (AP)			
S	60 mm	43 mm	65 mm	45 mm			
M	70 mm	48 mm </tr <tr> <td>L</td> <td>80 mm</td> <td>52 mm</td> <td>85 mm</td> <td>53 mm</td> </tr>	L	80 mm	52 mm	85 mm	53 mm
L	80 mm	52 mm	85 mm	53 mm			



# COPAL<sup>®</sup> knee moulds

## Filling volume

### FILLING VOLUME FOR THE FEMORAL COMPONENT

Size	Recommended number of packages of bone cement (1 x 40)
S	1
M	2
L	2

### FILLING VOLUME FOR THE TIBIAL COMPONENT

Size	Height	Recommended number of packages of bone cement (1 x 40)
S	Up to 32 mm	2
S	From 37 mm	3
M	Up to 27 mm	2
M	From 32 mm	3
L	Up to 22 mm	2
L	From 27 mm	3



[www.heraeus-medical.com](http://www.heraeus-medical.com)

66039571/  
06558